

## HDPE & PP Ergo Carboys

FOR GMP, PILOT PLANT & RESEARCH AND DEVELOPMENT  
Use with media, buffers, reagents, bulk API's and water



New! Space Saving Design, higher packing density of bulks

### Description

- HDPE and PP have broad chemical resistance, Also Available in Clear, Amber and Dark Amber.

Material HDPE, Polypropylene, Copolyester, Polycarb

- 20L, 40L and 75L Working Volume Capacity,
- Re-Usable with or w/o Disposable Liners
- FLEXI-CAP features connectors, dip tubes and QDC's.
- Ergonomic Easy to clean design
- Large (120mm) Diameter Mouth Lid

### CAUTION:

- Do not use for vacuum or pressure
- Do not use at temps higher than rating
- Do not use for storing oxidizing agents
- Do not flame carboy
- Do not use spigot as handle

### Ordering Options:



Dip tubes, Hose barbs, Quick Disconnect Fittings, sanitary flange lid, sterile or non-sterile as well by special order.

Autoclave	Polypropylene Only
Dimensions	See Sizing Charts Below
Temperature range	-4° C to + 50°C

### Applications

Reagents, Media, Buffers, Chromatography Solutions, High purity water, bulks, intermediates, long term storage. Shipping container




Clear Copolyester and Clear Polycarbonate (20L only)

Made in the USA



### Ordering Information:

Description - HDPE	Part Number	Lid Dim.	Pricing
20L Carboy , Versatile Cap	52- 3121-OEM	120mm	\$108.86
20L Carboy , Versatile Cap, w-spigot	52- 3221-OEM	120mm	\$168.26
40L Carboy , Versatile Cap	52- 5121-OEM	120mm	\$187.86
40L Carboy , Versatile Cap, w-spigot	52- 5221-OEM	120mm	\$243.56
75L Carboy , Versatile Cap, no spigot	52- 7121-OEM	120mm	\$361.56














Distributed by:  
**Optimum Processing Inc.**

For more information please go to: [www.marinscientific.com](http://www.marinscientific.com)

55 Mitchell Blvd. No. 23, Ca. 94903 Phone 415-461-7033 Fax 888-579-8824

To place an order email: [orders@opibioprocess.com](mailto:orders@opibioprocess.com)

**120mm cap adaptor options:**

 <p><b>205-2003-RLS</b> Cap, Closed Cap Size: 120mm Hose Barb Size: n/a No. of Quick Connects: n/a Qty/pk: 2</p>	 <p><b>205-4004-RLS</b> Adaptor, 1/2" HB x 3 Cap Size: 120mm Hose Barb Size: 1/2" No. of Quick Connects: 3 Qty/pk: 1</p>
 <p><b>205-2004-RLS</b> Cap, Open with closed adaptor Cap Size: 120mm Hose Barb Size: n/a No. of Quick Connects: n/a Qty/pk: 2</p>	 <p><b>205-4005-RLS</b> Adaptor, 3/4" HB x 3 Cap Size: 120mm Hose Barb Size: 3/4" No. of Quick Connects: 3 Qty/pk: 1</p>
 <p><b>205-4001-RLS</b> Adaptor, Closed Cap Size: 120mm Hose Barb Size: n/a No. of Quick Connects: n/a Qty/pk: 1</p>	 <p><b>205-4006-RLS</b> Adaptor, 1/2" HB x 4 Cap Size: 120mm Hose Barb Size: 1/2" No. of Quick Connects: 4 Qty/pk: 1</p>
 <p><b>205-4002-RLS</b> Adaptor, 1/2" HB x 2 Cap Size: 120mm Hose Barb Size: 1/2" No. of Quick Connects: 2 Qty/pk: 1</p>	 <p><b>205-4007-RLS</b> Adaptor, 3/4" HB x 4 Cap Size: 120mm Hose Barb Size: 3/4" No. of Quick Connects: 4 Qty/pk: 1</p>
 <p><b>205-4003-RLS</b> Adaptor, 3/4" HB x 2 Cap Size: 120mm Hose Barb Size: 3/4" No. of Quick Connects: 2 Qty/pk: 1</p>	 <p><b>205-4008-RLS</b> Adaptor, EZ Top® adaptor Cap Size: 120mm Hose Barb Size: n/a No. of Quick Connects: n/a Qty/pk: 1</p>
 <p><b>225-1001-RLS</b> Spigot, 2 pack Cap Size: n/a Hose Barb Size: n/a No. of Quick Connects: n/a Qty/pk: 2</p>	

HB = Hose Barb

**Lifting Safely Recommendations:**



**Distributed by:**  
**Optimum Processing Inc.**

For more information please go to: [www.marinscientific.com](http://www.marinscientific.com)

55 Mitchell Blvd. No. 23, Ca. 94903 Phone 415-461-7033 Fax 888-579-8824

To place an order email: [orders@opbioprocess.com](mailto:orders@opbioprocess.com)